



GLTP Rabbit pAb

Catalog No	YP-Ab-19241
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	GLTP
Immunogen	Synthesized peptide derived from human GLTP
Specificity	This antibody detects endogenous levels of GLTP at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Glycolipid transfer protein;GLTP;Glycolipid transfer protein (GLTP)
Observed Band	24 kDa
Calculated Molecular Weight	24 kDa
Cell Pathway	Cytoplasm .
Tissue Specificity	Detected in fibroblasts (at protein level). Detected in fibroblasts and in various cancer cell lines.
Function	Accelerates the intermembrane transfer of various glycolipids. Catalyzes the transfer of various glycosphingolipids between membranes but does not catalyze the transfer of phospholipids. May be involved in the intracellular translocation of glucosylceramides.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images