



IRF2BP1 Rabbit pAb

Catalog No	YP-Ab-19267
Isotype	IgG
Reactivity	Human,Mouse
Applications	WB
Gene Name	IRF2BP1
Immunogen	Synthesized peptide derived from human IRF2BP1
Specificity	This antibody detects endogenous levels of IRF2BP1 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	EC:2.3.2.27;Interferon regulatory factor 2-binding protein 1;IRF 2 binding protein 1;IRF 2BP1;IRF-2-binding protein 1;IRF2BP1;Interferon regulatory factor 2-binding protein 1 (IRF-2-binding protein 1) (IRF-2BP1) (Probable E3 ubiquitin-protein ligase IRF2BP1) (EC 6.3.2.-)
Observed Band	62-72 kDa
Calculated Molecular Weight	584 aa, 62 kDa
Cell Pathway	Nucleus .
Tissue Specificity	
Function	Acts as a transcriptional corepressor in a IRF2-dependent manner; this repression is not mediated by histone deacetylase activities. May act as an E3 ligase towards JDP2, enhancing its polyubiquitination. Represses ATF2-dependent transcriptional activation.
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images