



CSDE1 Rabbit pAb

Catalog No	YP-Ab-19312
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CSDE1 D1S155E KIAA0885 NRU UNR
Immunogen	Synthesized peptide derived from human CSDE1
Specificity	This antibody detects endogenous levels of CSDE1 at Human, Mouse,Rat
Formulation	
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	UNR;D1S155E;E1;CSDE1 D1S155E KIAA0885 NRU UNR;Cold shock domain-containing protein E1 (N-ras upstream gene protein) (Protein UNR)
Observed Band	89kDa
Calculated Molecular Weight	89kDa
Cell Pathway	Cytoplasm. Cytoplasm, Stress granule . Cytoplasm, P-body .
Tissue Specificity	
Function	RNA-binding protein involved in translationally coupled mRNA turnover . Implicated with other RNA-binding proteins in the cytoplasmic deadenylation/ translational and decay interplay of the FOS mRNA mediated by the major coding-region determinant of instability (mCRD) domain . Required for efficient formation of stress granules . ; (Microbial infection) Required for internal initiation of translation of human rhinovirus RNA.
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images