



UCHL5IP Rabbit pAb

Catalog No	YP-Ab-19323
Isotype	IgG
Reactivity	Human,Mouse
Applications	WB
Gene Name	HAUS7 UCHL5IP UIP1
Immunogen	Synthesized peptide derived from human UCHL5IP
Specificity	This antibody detects endogenous levels of UCHL5IP at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	UIP1;UCHL5IP;HAUS7 UCHL5IP UIP1;HAUS augmin-like complex subunit 7 (26S proteasome-associated UCH37-interacting protein 1) (UCHL5-interacting protein) (X-linked protein STS1769)
Observed Band	40kDa
Calculated Molecular Weight	40kDa
Cell Pathway	Cytoplasm, cytoskeleton, microtubule organizing center, centrosome . Cytoplasm, cytoskeleton, spindle . Localizes to interphase centrosomes and to mitotic spindle microtubules. .
Tissue Specificity	Detected in spleen, thymus, testis, ovary, small intestine and colon, with highest levels of expression in testis and ovary.
Function	Contributes to mitotic spindle assembly, maintenance of centrosome integrity and completion of cytokinesis as part of the HAUS augmin-like complex.
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images