



EML2 Rabbit pAb

Catalog No	YP-Ab-19360
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	EML2 EMAP2 EMAPL2
Immunogen	Synthesized peptide derived from human EML2. AA range:99-169
Specificity	This antibody detects endogenous levels of EML2 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from Rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	EML2 EMAP2 EMAPL2;Echinoderm microtubule-associated protein-like 2 (EMAP-2) (HuEMAP-2)
Observed Band	
Calculated Molecular Weight	71kD
Cell Pathway	Cytoplasm, cytoskeleton . Cytoplasm, cytoskeleton, spindle . Colocalizes with the microtubule cytoskeleton. Colocalizes with the mitotic spindle. .
Tissue Specificity	Ubiquitous.
Function	Tubulin binding protein that inhibits microtubule nucleation and growth, resulting in shorter microtubules.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images