



CRIP1 Rabbit pAb

Catalog No	YP-Ab-19363
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	CRIP1 CRIP CRP1
Immunogen	Synthesized peptide derived from human CRIP1. AA range:1-77
Specificity	This antibody detects endogenous levels of CRIP1 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CRIP1 CRIP CRP1;Cysteine-rich protein 1 (CRP-1) (Cysteine-rich heart protein) (CRHP) (hCRHP) (Cysteine-rich intestinal protein) (CRIP)
Observed Band	
Calculated Molecular Weight	8kD
Cell Pathway	
Tissue Specificity	
Function	Seems to have a role in zinc absorption and may function as an intracellular zinc transport protein.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images