



HDHC3 Rabbit pAb

Catalog No	YP-Ab-19390
Isotype	IgG
Reactivity	Human,Mouse
Applications	WB
Gene Name	HDHC3 MESH1
Immunogen	Synthesized peptide derived from human HDHC3. AA range:3-103
Specificity	This antibody detects endogenous levels of HDHC3 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	HDHC3 MESH1;Guanosine-3';5'-bis (diphosphate) 3'-pyrophosphohydrolase MESH1 (HD domain-containing protein 3) (Metazoan SpoT homolog 1) (MESH1) (Penta-phosphate guanosine-3'-pyrophosphohydrolase) (ppGpp) ase)
Observed Band	
Calculated Molecular Weight	20kD
Cell Pathway	
Tissue Specificity	
Function	ppGpp hydrolyzing enzyme involved in starvation response.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images