



RPL32 Rabbit pAb

Catalog No	YP-Ab-19394
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	RPL32 PP9932
Immunogen	Synthesized peptide derived from human RPL32. AA range:1-130
Specificity	This antibody detects endogenous levels of RPL32 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	RPL32 PP9932;60S ribosomal protein L32
Observed Band	
Calculated Molecular Weight	15kD
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images