



ZBTB24 Rabbit pAb

Catalog No	YP-Ab-19395
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	ZBTB24 KIAA0441 ZNF450
Immunogen	Synthesized peptide derived from human ZBTB24. AA range:429-499
Specificity	This antibody detects endogenous levels of ZBTB24 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ZBTB24 KIAA0441 ZNF450;Zinc finger and BTB domain-containing protein 24 (Zinc finger protein 450)
Observed Band	
Calculated Molecular Weight	77kD
Cell Pathway	Nucleus .
Tissue Specificity	Widely expressed , with highest levels in naive B-cells.
Function	May be involved in BMP2-induced transcription.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images