



GYS2 rabbit pAb

Catalog No	YP-Ab-10024
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	GYS2
Immunogen	Synthesized peptide derived from human GYS2 AA range: 26-76
Specificity	This antibody detects endogenous levels of human GYS2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	glycogen synthase 2 (liver);Glycogen synthase 2;EC:2.4.1.11;GYS2;Glycogen [starch] synthase;liver (EC 2.4.1.11)
Observed Band	81 kDa
Calculated Molecular Weight	703 aa, 81 kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
