



PGTB2 rabbit pAb

Catalog No	YP-Ab-10031
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	RABGGTB GGTB
Immunogen	Synthesized peptide derived from human PGTB2 AA range: 64-114
Specificity	This antibody detects endogenous levels of human PGTB2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PGTB2;RABGGTB GGTB;Geranylgeranyl transferase type-2 subunit beta (EC 2.5.1.60) (Geranylgeranyl transferase type II subunit beta) (GGTase-II-beta) (Rab geranyl-geranyltransferase subunit beta) (Rab GG transferase beta)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
