



ZF64B rabbit pAb

Catalog No	YP-Ab-10038
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	ZFP64 ZNF338
Immunogen	Synthesized peptide derived from human ZF64B AA range: 292-342
Specificity	This antibody detects endogenous levels of human ZF64B
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ZF64B;ZFP64 ZNF338;Zinc finger protein 64 homolog;isoforms 3 and 4 (Zfp-64) (Zinc finger protein 338)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
