



EI2BA rabbit pAb

Catalog No	YP-Ab-10043
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	EIF2B1 EIF2BA
Immunogen	Synthesized peptide derived from human EI2BA AA range: 53-103
Specificity	This antibody detects endogenous levels of human EI2BA
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	EI2BA_HUMAN;eIF-2B GDP-GTP exchange factor subunit alpha;EIF2B;Eif2b1; EIF2BA;Eukaryotic translation initiation factor 2B subunit 1 alpha 26kDa; Eukaryotic translation initiation factor 2B subunit alpha;Translation initiation factor eIF-2B subunit alpha;EI2BA;EIF2B1 EIF2BA;Translation initiation factor eIF-2B subunit alpha (eIF-2B GDP-GTP exchange factor subunit alpha)
Observed Band	34kD
Calculated Molecular Weight	34 kDa;
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

