



RBM14 rabbit pAb

Catalog No	YP-Ab-10076
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	RBM14 SIP
Immunogen	Synthesized peptide derived from human RBM14 AA range: 465-515
Specificity	This antibody detects endogenous levels of human RBM14
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	RNA-binding protein 14;RNA-binding motif protein 14;RNA binding protein 14; PSP2;Paraspeckle protein 2;RBM14;RBM14 SIP;RNA-binding protein 14 (Paraspeckle protein 2) (PSP2) (RNA-binding motif protein 14) (RRM-containing coactivator activator/modulator) (Synaptotagmin-interacting protein) (SYT-interacting protein)
Observed Band	70-75 kDa
Calculated Molecular Weight	70 kDa
Cell Pathway	Nucleus . Nucleus, nucleolus . Cytoplasm . In punctate subnuclear structures often located adjacent to splicing speckles, called paraspeckles (PubMed:11790299). Cytoplasmic localization is crucial for its function in suppressing the formation of aberrant centriolar protein complexes (PubMed:25385835). .
Tissue Specificity	Expressed in all tissues tested, including brain, heart, skeletal muscle, colon, thymus, spleen, kidney, liver, small intestine, placenta, lung and peripheral blood lymphocytes.
Function	
Background	

**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images