



AMFR rabbit pAb

Catalog No	YP-Ab-10091
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	AMFR RNF45
Immunogen	Synthesized peptide derived from human AMFR AA range: 572-622
Specificity	This antibody detects endogenous levels of human AMFR
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	AMFR;E3 ubiquitin-protein ligase AMFR;EC:2.3.2.36;GP78;RNF 45;AMFR RNF45;E3 ubiquitin-protein ligase AMFR (EC 6.3.2.-) (Autocrine motility factor receptor;isoform 2) (AMF receptor;isoform 2) (RING finger protein 45) (gp78)
Observed Band	73 kDa
Calculated Molecular Weight	643 aa, 73 kDa
Cell Pathway	Endoplasmic reticulum membrane ; Multi-pass membrane protein . Palmitoylation promotes localization to the peripheral endoplasmic reticulum. .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
