



CC019 rabbit pAb

Catalog No	YP-Ab-10096
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	C3orf19
Immunogen	Synthesized peptide derived from human CC019 AA range: 14-64
Specificity	This antibody detects endogenous levels of human CC019
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CC019;C3orf19;Uncharacterized protein C3orf19
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus . In the nucleus of neuroblastoma cells, it colocalizes with EIF4A3, a component of exon junction complex. .
Tissue Specificity	Widely expressed.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
