



ME3L2 rabbit pAb

Catalog No	YP-Ab-10109
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	MEIS3P2
Immunogen	Synthesized peptide derived from human ME3L2 AA range: 200-250
Specificity	This antibody detects endogenous levels of human ME3L2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ME3L2;MEIS3P2;Putative homeobox protein Meis3-like 2
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
