



KDM7 rabbit pAb

Catalog No	YP-Ab-10116
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	JHDM1D KDM7 KIAA1718
Immunogen	Synthesized peptide derived from human KDM7 AA range: 217-267
Specificity	This antibody detects endogenous levels of human KDM7
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	KDM7;JHDM1D KDM7 KIAA1718;Lysine-specific demethylase 7 (EC 1.14.11.-) (JmjC domain-containing histone demethylation protein 1D)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus.
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
