



SERF1 rabbit pAb

Catalog No	YP-Ab-10121
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	SERF1A FAM2A SERF1 SMAM1; SERF1B FAM2B SERF1 SMAM1
Immunogen	Synthesized peptide derived from human SERF1 AA range: 10-60
Specificity	This antibody detects endogenous levels of human SERF1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SERF1;SERF1A FAM2A SERF1 SMAM1;SERF1B FAM2B SERF1 SMAM1; Small EDRK-rich factor 1 (Protein 4F5) (h4F5) (SMA modifier 1)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm, cytosol . Nucleus .
Tissue Specificity	Isoform Long is predominantly expressed in heart, brain and skeletal muscle. Isoform Short and Isoform Long are expressed throughout the central nervous system, including spinal cord.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
