



PDE6C rabbit pAb

Catalog No	YP-Ab-10122
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	PDE6C PDEA2
Immunogen	Synthesized peptide derived from human PDE6C AA range: 474-524
Specificity	This antibody detects endogenous levels of human PDE6C
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	5"-cyclic phosphodiesterase subunit alpha";cGMP phosphodiesterase 6C;COD4; Cone cGMP specific 3' 5' cyclic phosphodiesterase subunit alpha';Cone cGMP-specific 3";PDE 6C;PDE6 alpha prime;PDE6 alpha';PDE6C;PDE6C_HUMAN; PDEA2;Phosphodiesterase 6C cGMP specific cone alpha prime;PDE6C PDEA2; Cone cGMP-specific 3';5'-cyclic phosphodiesterase subunit alpha' (EC 3.1.4.35) (cGMP phosphodiesterase 6C)
Observed Band	99kD
Calculated Molecular Weight	99kD
Cell Pathway	Cell membrane ; Lipid-anchor ; Cytoplasmic side .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
