



NARG2 rabbit pAb

Catalog No	YP-Ab-10125
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	NARG2 BRCC1 UNQ3101/PRO10100
Immunogen	Synthesized peptide derived from human NARG2 AA range: 474-524
Specificity	This antibody detects endogenous levels of human NARG2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	BRCC1;ICE2;Interactor of little elongator complex ELL subunit 2;Little elongation complex subunit 2;NMDA receptor-regulated protein 2;NARG2; NARG2 BRCC1 UNQ3101/PRO10100
Observed Band	140 kDa
Calculated Molecular Weight	982 aa, 109 kDa
Cell Pathway	Nucleus . Colocalizes with COIL in subnuclear Cajal and histone locus bodies. Translocates in the LEC complex to Cajal and histone locus bodies at snRNA genes in a ICE1-dependent manner. Associates to transcriptionally active chromatin at snRNA genes.
Tissue Specificity	Expressed at low levels in lung and testis.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images