



CK037 rabbit pAb

Catalog No	YP-Ab-10129
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	LINC00167 C11orf37 NCRNA00167
Immunogen	Synthesized peptide derived from human CK037 AA range: 55-105
Specificity	This antibody detects endogenous levels of human CK037
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CK037;LINC00167 C11orf37 NCRNA00167;Putative uncharacterized protein encoded by LINC00167
Observed Band	16kD
Calculated Molecular Weight	16kD
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

