



TM74B rabbit pAb

Catalog No	YP-Ab-10132
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	TMEM74B C20orf46
Immunogen	Synthesized peptide derived from human TM74B AA range: 52-102
Specificity	This antibody detects endogenous levels of human TM74B
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TM74B;TMEM74B C20orf46;Transmembrane protein 74B
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Membrane ; Multi-pass membrane protein .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
