



MMSA rabbit pAb

Catalog No	YP-Ab-10162
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ALDH6A1 MMSDH
Immunogen	Synthesized peptide derived from human MMSA AA range: 220-270
Specificity	This antibody detects endogenous levels of human MMSA
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Antigen MMSA 1;Asp His His Cys domain containing protein 9;CGI89;CXorf11; DHHC-9;DHHC9;MMSA1;MRXSZ;Palmitoyltransferase ZDHHC9; ZDHC9_HUMAN;ZDHHC 9;ZDHHC10;ZDHHC9;Zinc finger DHHC domain containing 10;Zinc finger DHHC domain containing protein 9;Zinc finger DHHC domain-containing protein 9;Zinc finger DHHC type containing 9;Zinc finger protein 379;Zinc finger protein 380;ZNF379;ZNF380;MMSA;ALDH6A1 MMSDH; Methylmalonate-semialdehyde dehydrogenase [acylating];mitochondrial (MMSDH) (Malonate-semialdehyde dehydrogenase [acylating]) (EC 1.2.1.18) (EC 1.2.1.27) (Aldehyde dehydrogenase family 6 member A1)
Observed Band	41kD
Calculated Molecular Weight	41kDa
Cell Pathway	Mitochondrion.
Tissue Specificity	
Function	
Background	

**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
