



ZNF2 rabbit pAb

Catalog No	YP-Ab-10170
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	ZNF2 ZNF661
Immunogen	Synthesized peptide derived from human ZNF2 AA range: 25-75
Specificity	This antibody detects endogenous levels of human ZNF2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	A1 5;Zfp661;zinc finger protein 2;Zinc finger protein 2.2;Zinc finger protein 661; ZNF2;ZNF661;ZNF2 ZNF661;Zinc finger protein 2 (Zinc finger protein 2.2) (Zinc finger protein 661)
Observed Band	49 kDa
Calculated Molecular Weight	49 kDa
Cell Pathway	Nucleus .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
