



PRR11 rabbit pAb

Catalog No	YP-Ab-10190
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	PRR11
Immunogen	Synthesized peptide derived from human PRR11 AA range: 282-332
Specificity	This antibody detects endogenous levels of human PRR11
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Proline rich protein 11;proline rich 11;PRR11;Proline-rich protein 11
Observed Band	38 kDa
Calculated Molecular Weight	360 aa, 40 kDa
Cell Pathway	Cytoplasm . Nucleus .
Tissue Specificity	Ubiquitously expressed.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

