



PGAM1 rabbit pAb

Catalog No	YP-Ab-10237
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PGAM1 PGAMA CDABP0006
Immunogen	Synthesized peptide derived from human PGAM1 AA range: 14-64
Specificity	This antibody detects endogenous levels of human PGAM1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PGAM-B;PGAMA;PGAM 1;EC:5.4.2.11;CDABP0006;PGAM1;PGAM1 PGAMA CDABP0006;Phosphoglycerate mutase 1 (EC 3.1.3.13) (EC 5.4.2.1) (EC 5.4.2.4) (BPG-dependent PGAM 1) (Phosphoglycerate mutase isozyme B) (PGAM-B)
Observed Band	29 kDa
Calculated Molecular Weight	254 aa, 29 kDa
Cell Pathway	
Tissue Specificity	Expressed in the liver and brain. Not found in the muscle.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
