



KIF3B rabbit pAb

Catalog No	YP-Ab-10238
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	KIF3B KIAA0359
Immunogen	Synthesized peptide derived from human KIF3B AA range: 256-306
Specificity	This antibody detects endogenous levels of human KIF3B
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	KIF3B;KIF3B KIAA0359;Kinesin-like protein KIF3B (HH0048) (Microtubule plus end-directed kinesin motor 3B)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm, cytoskeleton . Cell projection, cilium . Cell projection, dendritic spine .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
