



PHKG2 rabbit pAb

Catalog No	YP-Ab-10245
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PHKG2
Immunogen	Synthesized peptide derived from human PHKG2 AA range: 166-216
Specificity	This antibody detects endogenous levels of human PHKG2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	EC:2.7.11.19;GSD9C;PHK gamma T;PHK-gamma-LT;PHK-gamma-T;PHKG2; Phosphorylase b kinase gamma catalytic chain;testis/liver isoform (PHK-gamma-T) (EC 2.7.11.19) (PSK-C3) (Phosphorylase kinase subunit gamma-2) (Serine/ threonine-protein kinase PHKG2) (EC 2.7.11.1) (E
Observed Band	46 kDa
Calculated Molecular Weight	46 kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
