



GRD2I rabbit pAb

Catalog No	YP-Ab-10254
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	GRID2IP
Immunogen	Synthesized peptide derived from human GRD2I AA range: 43-93
Specificity	This antibody detects endogenous levels of human GRD2I
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	GRD2I;GRID2IP;Delphilin (Glutamate receptor;ionotropic;delta 2-interacting protein 1)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cell junction, synapse, postsynaptic cell membrane .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
