



PHF2 rabbit pAb

Catalog No	YP-Ab-10260
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	PHF2 KIAA0662
Immunogen	Synthesized peptide derived from human PHF2 AA range: 74-124
Specificity	This antibody detects endogenous levels of human PHF2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	GRC5;JHDM1E;KIAA0662;PHD finger protein 2;PHF2;PHF2 KIAA0662;Lysine-specific demethylase PHF2 (EC 1.14.11.-) (GRC5) (PHD finger protein 2)
Observed Band	150 kDa
Calculated Molecular Weight	1101 aa, 121 kDa
Cell Pathway	Nucleus, nucleolus . Chromosome, centromere, kinetochore .
Tissue Specificity	Widely expressed, including in liver (at protein level).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
