



# ENTP3 rabbit pAb

|                                    |                                                                                                                                                                                                                         |
|------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Catalog No</b>                  | YP-Ab-10265                                                                                                                                                                                                             |
| <b>Isotype</b>                     | IgG                                                                                                                                                                                                                     |
| <b>Reactivity</b>                  | Human;Rat;Mouse;                                                                                                                                                                                                        |
| <b>Applications</b>                | WB                                                                                                                                                                                                                      |
| <b>Gene Name</b>                   | ENTPD3 CD39L3                                                                                                                                                                                                           |
| <b>Immunogen</b>                   | Synthesized peptide derived from human ENTP3 AA range: 298-348                                                                                                                                                          |
| <b>Specificity</b>                 | This antibody detects endogenous levels of human ENTP3                                                                                                                                                                  |
| <b>Formulation</b>                 | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                                                                                                                                 |
| <b>Source</b>                      | Polyclonal, Rabbit,IgG                                                                                                                                                                                                  |
| <b>Purification</b>                | The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.                                                                                                               |
| <b>Dilution</b>                    | WB 1:1000-2000                                                                                                                                                                                                          |
| <b>Concentration</b>               | 1 mg/ml                                                                                                                                                                                                                 |
| <b>Purity</b>                      | ≥90%                                                                                                                                                                                                                    |
| <b>Storage Stability</b>           | -20°C/1 year                                                                                                                                                                                                            |
| <b>Synonyms</b>                    | ENTP3;ENTPD3 CD39L3;Ectonucleoside triphosphate diphosphohydrolase 3 (NTPDase 3) (EC 3.6.1.5) (CD39 antigen-like 3) (Ecto-ATP diphosphohydrolase 3) (Ecto-ATPDase 3) (Ecto-ATPase 3) (Ecto-apyrase 3) (HB6)             |
| <b>Observed Band</b>               |                                                                                                                                                                                                                         |
| <b>Calculated Molecular Weight</b> |                                                                                                                                                                                                                         |
| <b>Cell Pathway</b>                | Cell membrane ; Multi-pass membrane protein .                                                                                                                                                                           |
| <b>Tissue Specificity</b>          | Expressed in adult brain, pancreas, spleen and prostate (PubMed:9676430). Moderate or low expression is seen in most tissues (PubMed:9676430). Not expressed in liver and peripheral blood leukocytes (PubMed:9676430). |
| <b>Function</b>                    |                                                                                                                                                                                                                         |
| <b>Background</b>                  |                                                                                                                                                                                                                         |
| <b>matters needing attention</b>   | Avoid repeated freezing and thawing!                                                                                                                                                                                    |
| <b>Usage suggestions</b>           | This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.                                                                                       |



## Products Images