




mStrawberry Polyclonal Antibody

Catalog No	YP-Ab-10275
Isotype	IgG
Reactivity	Species independent
Applications	WB
Gene Name	
Immunogen	Recombinant Protein of mStrawberry
Specificity	The antibody detects mStrawberry & mStrawberry Tag protein.
Formulation	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.
Dilution	WB: 1:2000-5000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	mStrawberry;YM3234
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	mStrawberry is a bright, monomeric red fluorescent protein. mStrawberry excitation and emission maxima are 574 and 596 nm, respectively.
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of mStawberry recombinant protein, diluted at 1:5000. Secondary antibody(catalog#:RS0002) was diluted at 1:20000