




GST-Tag Polyclonal Antibody

Catalog No	YP-Ab-10280
Isotype	IgG
Reactivity	Species independent
Applications	WB;ELISA
Gene Name	
Immunogen	Recombinant Protein of GST-Tag
Specificity	The antibody detects C-terminal, internal, and N-terminal GST fusion proteins.
Formulation	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.
Dilution	WB: 1:5000-10000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	GST-Tag;YM3205
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	GST (Glutathione S-Transferase) is a widely used protein tag encoded by the Schistosoma japonicum. GST provides both an easily detectable tag and a simple purification process with little effect on the biological function of the protein of interest. Antibodies to GST are useful for checking protein expression both in plaques and on Western blots as well as for immunoaffinity purification of proteins expressed in cells.
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

0.5ug GST fusion protein+ Primary antibody dilution at
1) 1:5000 2) 1:20000. Secondary
antibody(catalog#:RS0002) was diluted at 1:20000