



BFP Polyclonal Antibody

Catalog No	YP-Ab-10288
Isotype	IgG
Reactivity	Species independent
Applications	WB
Gene Name	
Immunogen	Recombinant Protein of BFP
Specificity	The antibody detects EBFP & EBFP Tag protein.
Formulation	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.
Dilution	WB: 1:5000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	BFP;YM3238
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	Blue Fluorescent Protein (BFP) emits a strong blue fluorescence. EBFP can be used as protein "tags" to study the sub-cellular localization of proteins and/or their translocation upon stimulation or as markers for transfection in transient and stable expression systems.
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of EBFP recombinant protein, diluted at 1:5000. Secondary antibody(catalog#:RS0002) was diluted at 1:20000