




Actin beta rabbit pAb(Caenorhabditis elegans)

Catalog No	YP-Ab-10297
Isotype	IgG
Reactivity	C.elegans
Applications	WB
Gene Name	
Immunogen	Purified recombinant Full length Actin beta (Caenorhabditis elegans)
Specificity	This antibody detects endogenous levels of Actin beta (Caenorhabditis elegans)
Formulation	Liquid in PBS containing 50% glycerol and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Actin beta(Caenorhabditis elegans);Actin beta;Beta Actin;βActin
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of C.elegans cells using Actin beta rabbit pAb(Caenorhabditis elegans). Secondary antibody(catalog#:RS0002) was diluted at 1:20000