




Tubulin beta rabbit pAb(Caenorhabditis elegans)

Catalog No	YP-Ab-10300
Isotype	IgG
Reactivity	C.elegans
Applications	WB
Gene Name	
Immunogen	Purified recombinant Full lenth Tubulin beta (Caenorhabditis elegans)
Specificity	This antibody detects endogenous levels ofTubulin beta (Caenorhabditis elegans)
Formulation	Liquid in PBS containing 50% glycerol and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Tubulin beta 3 chain;Tubulin beta III;β III tubulin;βIII tubulin;TUBB4;Tubulin beta; Beta Tubulin
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of C.elegans cells using Tubulin beta rabbit pAb(Caenorhabditis elegans). Secondary antibody(catalog#:RS0002) was diluted at 1:20000