




GAPDH rabbit pAb(Drosophilidae)

Catalog No	YP-Ab-10302
Isotype	IgG
Reactivity	Drosophilidae
Applications	WB
Gene Name	
Immunogen	Purified recombinant Full lenth GAPDH (Drosophiladaes)
Specificity	This antibody detects endogenous levels of GAPDH (Drosophiladaes)
Formulation	Liquid in PBS containing 50% glycerol and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	GAPD;Glyceraldehyde-3-phosphate dehydrogenase;GAPDH
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of Drosophilidae cells using GAPDH rabbit pAb(Drosophilidae). Secondary antibody(catalog#:RS0002) was diluted at 1:20000