




Tubulin alpha rabbit pAb(Yeast)

Catalog No	YP-Ab-10304
Isotype	IgG
Reactivity	Yeast
Applications	WB
Gene Name	
Immunogen	Purified recombinant Full lenth Tubulin alpha (Yeast)
Specificity	This antibody detects endogenous levels of Tubulin alpha (Yeast)
Formulation	Liquid in PBS containing 50% glycerol and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TUBA1B;Tubulin alpha 1B chain;Alpha-tubulin ubiquitous;EC:3.6.5.-;Tubulin alpha
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of Yeast cells using Tubulin alpha rabbit pAb(Yeast). Secondary antibody(catalog#:RS0002) was diluted at 1:20000