



SARS-COV-2 Spike Protein S1 Rabbit pAb

Catalog No	YP-Ab-10310
Isotype	IgG
Reactivity	Human virus
Applications	WB
Gene Name	SARS-COV-2 Spike Protein S1
Immunogen	Synthesized peptide derived from SARS-COV-2 Spike Protein S1
Specificity	This antibody detects endogenous levels of SARS-COV-2 Spike Protein S1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.17% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SARS-CoV-2 Spike Glycoprotein (428-56 aa);SARS-COV-2 Spike;SARS-COV-2 Spike Protein S1;YM33057
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	The spike (S) glycoprotein of coronaviruses is known to be essential in the binding of the virus to the host cell at the advent of the infection process. Most notable is severe acute respiratory syndrome (SARS). The severe acute respiratory syndrome-coronavirus (SARS-CoV) spike (S) glycoprotein alone can mediate the membrane fusion required for virus entry and cell fusion. It is also a major immunogen and a target for entry inhibitors. It's been reported that 2019-nCoV can infect the human respiratory epithelial cells through interaction with the human ACE2 receptor.


**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western blot analysis of recombinant SARS-CoV-2 Spike Protein using Rabbit pAb diluted at 1:10,000.