




mScarlet Rabbit pAb


Catalog No	YP-Ab-10316
Isotype	IgG
Reactivity	All
Applications	IHC;WB
Gene Name	
Immunogen	Synthesized peptide derived from mScarlet
Specificity	This antibody detects endogenous levels of mScarlet
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.80% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	IHC-p1:50-200 ,WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	mScarlet;YM33120
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	mScarlet: a bright monomeric red fluorescent protein for cellular imaging.
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of mScarlet recombinant transfected 293T Cell Lysate using mScarlet Rabbit Polyclonal Ab diluted at 1:2,000.



Immunohistochemical analysis of paraffin-embedded human tonsil. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA, pH9.0 was used for antigen retrieval. 3, Secondary antibody was diluted at 1:200(room temperature, 30min).