



ARHGAP42 (phospho-Tyr376) rabbit pAb

Catalog No	YP-Ab-10358
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	ARHGAP42
Immunogen	Synthesized phosho peptide around human ARHGAP42 (Tyr376)
Specificity	This antibody detects endogenous levels of Human Mouse ARHGAP42 (phospho-Tyr376)
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ARHGAP42;ARHGAP42 (Tyr376);Rho GTPase-activating protein 42 (Rho GTPase-activating protein 10-like) (Rho-type GTPase-activating protein 42)
Observed Band	95kD
Calculated Molecular Weight	95kD
Cell Pathway	
Tissue Specificity	Highly and selectively expressed in smooth muscle cells.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
