



BACH1/BRIP1 (phospho-Thr1133) rabbit pAb

Catalog No	YP-Ab-10361
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	BRIP1 BACH1 FANCI
Immunogen	Synthesized phospho peptide around human BACH1 and BRIP1 (Thr1133)
Specificity	This antibody detects endogenous levels of Human BACH1 or BRIP1 (phospho-Thr1133)
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	BACH 1;BACH-1;BACH1/BRIP1;BRIP1 BACH1 FANCI;BACH1/BRIP1 (Thr1133);Fanconi anemia group J protein (Protein FANCI) (EC 3.6.4.13) (ATP-dependent RNA helicase BRIP1) (BRCA1-associated C-terminal helicase 1) (BRCA1-interacting protein C-terminal helicase 1) (BRCA1-interacting protein 1)
Observed Band	100-110 kDa
Calculated Molecular Weight	82 kDa
Cell Pathway	Nucleus . Cytoplasm .
Tissue Specificity	Ubiquitously expressed, with highest levels in testis.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images