



BCL9L (phospho-Ser915) rabbit pAb

Catalog No	YP-Ab-10363
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	BCL9L DLNB11
Immunogen	Synthesized phosho peptide around human BCL9L (Ser915)
Specificity	This antibody detects endogenous levels of Human BCL9L (phospho-Ser915)
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Nuclear Co-Factor Of Beta-Catenin Signalling;DLNB11;BCL9-like protein;B-cell lymphoma 9-like protein;B-cell CLL/lymphoma 9-like protein;BCL9L;BCL9L DLNB11;BCL9L (Ser915);B-cell CLL/lymphoma 9-like protein (B-cell lymphoma 9-like protein) (BCL9-like protein) (Protein BCL9-2)
Observed Band	220 kDa
Calculated Molecular Weight	157 kDa
Cell Pathway	Nucleus .
Tissue Specificity	Expressed in breast, ductal and invasive ductal carcinomas of the breast, sporadic colorectal adenomas and carcinomas (at protein level). Expressed in fetal brain. Expressed in lung, amygdala, eye, prostate, pancreatic and prostate cancers, head and neck tumors and embryonal tumor.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images