



CTDSPL2 (phospho-Ser104) rabbit pAb

Catalog No	YP-Ab-10364
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	CTDSPL2 HSPC058 HSPC129
Immunogen	Synthesized phosho peptide around human CTDSPL2 (Ser104)
Specificity	This antibody detects endogenous levels of Human CTDSPL2 (phospho-Ser104)
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CTD small phosphatase-like protein 2;CTDSP like 2;CTDSP-like 2;EC:3.1.3.-; HSPC058;CTDSPL2;CTDSPL2 HSPC058 HSPC129;CTDSPL2 (Ser104);CTD small phosphatase-like protein 2 (CTDSP-like 2) (EC 3.1.3.-)
Observed Band	67 kDa
Calculated Molecular Weight	466 aa, 53 kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
