



DNAJC2/MPP11 (phospho-Ser47) rabbit pAb

Catalog No	YP-Ab-10367
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	DNAJC2 MPHOSPH11 MPP11 ZRF1
Immunogen	Synthesized phosho peptide around human DNAJC2 and MPP11 (Ser47)
Specificity	This antibody detects endogenous levels of Human Mouse Rat DNAJC2/MPP11 (phospho-Ser47)
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	DNJC2;DNAJC2 MPHOSPH11 MPP11 ZRF1;DNAJC2/MPP11 (Ser47);DnaJ homolog subfamily C member 2 (M-phase phosphoprotein 11) (Zuotin-related factor 1)
Observed Band	68kD
Calculated Molecular Weight	68kD
Cell Pathway	Nucleus . Cytoplasm, cytosol .
Tissue Specificity	Widely expressed.
Function	PTM:Phosphorylated in M (mitotic) phase.,similarity:Contains 1 J domain., similarity:Contains 2 SANT domains.,
Background	This gene is a member of the M-phase phosphoprotein (MPP) family. The gene encodes a phosphoprotein with a J domain and a Myb DNA-binding domain which localizes to both the nucleus and the cytosol. The protein is capable of forming a heterodimeric complex that associates with ribosomes, acting as a molecular chaperone for nascent polypeptide chains as they exit the ribosome. This protein was identified as a leukemia-associated antigen and expression of the gene is upregulated in leukemic blasts. Also, chromosomal aberrations



involving this gene are associated with primary head and neck squamous cell tumors. This gene has a pseudogene on chromosome 6. Alternatively spliced variants which encode different protein isoforms have been described. [provided by RefSeq, Jul 2008],

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images