




MFF (phospho-Ser146) rabbit pAb


Catalog No	YP-Ab-10380
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB;ELISA;IHC
Gene Name	
Immunogen	Synthesized phosho peptide around human MFF (Ser146)
Specificity	This antibody detects endogenous levels of Human Mouse Rat MFF (phospho-Ser146)
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000;IHC-p 1:50-300; ELISA 2000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	AD033;C2orf33;GL004;Mitochondrial fission factor;MFF;MFF (Ser146)
Observed Band	37kD
Calculated Molecular Weight	37kD
Cell Pathway	Mitochondrion outer membrane ; Single-pass type IV membrane protein . Peroxisome . Cytoplasmic vesicle, secretory vesicle, synaptic vesicle .
Tissue Specificity	Highly expressed in heart, kidney, liver, brain, muscle, and stomach.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from U2OS cells, primary antibody was diluted at 1:1000, 4° over night



Immunohistochemical analysis of paraffin-embedded human liver cancer. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA, pH9.0 was used for antigen retrieval. 3, Secondary antibody was diluted at 1:200(room temperature, 45min).