



RSK2 (phospho-Ser227) rabbit pAb

Catalog No	YP-Ab-10403
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	
Immunogen	Synthesized phosho peptide around human RSK2 (Ser227)
Specificity	This antibody detects endogenous levels of Human Mouse Rat RSK2 (phospho-Ser227)
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	p90RSK3 (phospho Ser227);Phospho RSK2 (Ser227);Phospho-RPS6KA3; Phospho-RSK2;Phospho-RSK2 (Ser227);RSK2;RSK2 (Ser227)
Observed Band	84 kDa
Calculated Molecular Weight	740 aa, 84 kDa
Cell Pathway	Nucleus. Cytoplasm.
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
