



Semaphorin 4B (phospho-Ser825) rabbit pAb

Catalog No	YP-Ab-10406
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	
Immunogen	Synthesized phosho peptide around human Semaphorin 4B (Ser825)
Specificity	This antibody detects endogenous levels of Human Mouse Rat Semaphorin 4B (phospho-Ser825)
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Semaphorin 4B;Semaphorin 4B (Ser825)
Observed Band	95kD
Calculated Molecular Weight	95kD
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images