



Desmoyokin (Phospho Ser5782) rabbit pAb

Catalog No	YP-Ab-10421
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB; ELISA
Gene Name	AHNAK PM227
Immunogen	Synthesized peptide derived from human Desmoyokin (Phospho Ser5782)
Specificity	This antibody detects endogenous levels of Human Desmoyokin (Phospho Ser5782)
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Desmoyokin;AHNAK PM227;Desmoyokin (Phospho Ser5782);Neuroblast differentiation-associated protein AHNAK (Desmoyokin)
Observed Band	647kD
Calculated Molecular Weight	647kD
Cell Pathway	Nucleus.
Tissue Specificity	
Function	function:May be required for neuronal cell differentiation.,similarity:Contains 1 PDZ (DHR) domain.,subunit:Interacts with DYSF; the interaction is direct and Ca(2+)-independent.,
Background	function:May be required for neuronal cell differentiation.,similarity:Contains 1 PDZ (DHR) domain.,subunit:Interacts with DYSF; the interaction is direct and Ca(2+)-independent.,
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images